

Soil respiration measurement using Micro-Oxymax Respirometer

Basal respiration

- I. Soil preparation
 - a. 50g (dry weight) soil in 250ml flask¹ (with GL 45 fitting)
 - b. Setup the (water) bottles and tubing to Extension Interfaces
 - c. Weight wet soil based on moisture content measurement into numbered (according to channel number) Duran bottles.
 - d. Assemble as bottles are weighed except a tubing to Extension Interface to equilibrate with atmosphere
 - e. After all bottles are weighed, finish setup making sure O-rings and metal plates connects are intact.
 - f. Pre-incubate for 12-24 hours at incubation temperature (RT)²
- II. Setting up the Micro Oxymax Respirometer
 - a. Calibration (Tools-Calibration)
 - i. Offsetting gases other than O₂ (bottle 1 – lime soda column)
 - ii. Offsetting O₂ (bottle 2 – calibration gas cylinder³)
 - iii. Gaining O₂ (bottle 1)
 - iv. Gaining gases other than O₂ (bottle 2⁴)
 - b. Leakage check (Tools-Utilities)
 - i. Click on Leakage button for automatic testing (pp. 76 for source of leakage)
 - c. Setup the experiment (Experiment-Setup)
 - i. Can be done while leakage test is on
 - ii. Setup tab
 1. **Channels** – End channel – # channels to use
 2. Mark ‘Auto volume Measurement’, ‘O₂ Consumption Positive’
 3. **Timing** – Sample Interval(Hours) – 2.5hours
 4. **Data Units** – Normalization Unit – g
 - iii. Chamber Setup tab
 1. Type ‘Normalization Units’ and ‘Channel Label’
 2. ‘Volume’ and ‘Leakage’ will be determined by test
 - e. Start the incubation by clicking Run

¹ If available soil samples are not enough, try with as much as possible but be consistent with all samples within an experiment.

² Pre-incubation is needed for restoration of metabolic equilibrium of the population after storing at 4C (12hr) and re-wetting (24hr).

³ Open the cylinder first, open regulator knob until hearing hissing sound, then adjust pressure to 5psi.

⁴ Current composition (H₂S 179.8ppm, H₂ 1801ppm, CO₂ 0.906%, CO 0.9015%, CH₄ 0.9054%) – balanced by N₂ (cylinder # SG9147282, reference # 83-124147696-1). It should be modified in System Properties at Sensors tab (pp. 70)

- III. Incubation condition & measurement
 - a. 23C (or RT), 2 ½ hour interval⁵ for 120 hours⁶
 - b. Humidifying bottles (50ml with GL32 cap) for each sample to maintain moisture level in gas after going through drier (Figure 1)
 - c. Collect final results (accumulation) and graph in both rate ($\mu\text{l O}_2 \text{ g}^{-1} \text{ soil h}$) and accumulation ($\mu\text{l O}_2 \text{ g}^{-1} \text{ soil}$) for further inference
- IV. After incubation
 - a. Disconnect Duran bottles, empty soil, wash, oven-dry and cap with aluminum foil for long term storage; otherwise, re-use them for next incubation.
 - b. Leave water bottles connected for immediate next incubation; otherwise cover with the original orange caps

⁵ Or shortest intervals possible. (# channel + 1) * 6 min

⁶ 72 hours by Chuck, 150 hours by Hollender *et al* 2003